

Product Datasheet

Anti-Norovirus (GI)
V4925

Product Name	Anti-Norovirus (GI)
Catalog Number	V4925
Tested Applications	LF, ELISA, WB
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with GII
Specificity	GI capsid protein
Matched Pairs	LF IA pairs (capture-conjugate) : V4925-V4922
Shelf Life	Four years from date of manufacture