

Product Datasheet

Anti-Norovirus (GII) V4933

Product Name	Anti-Norovirus (GII)
Catalog Number	V4933
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with GI
Specificity	GII capsid protein
Matched Pairs	ELISA pairs (capture-conjugate): V4932-V4933 Lateral Flow pairs (capture-conjugate): V4935-V4933, V4935-V4935
Shelf Life	4 years