

Product Datasheet

Anti-Human CD13 Monoclonal Antibody, Azide Free Clone B-F10 854.100.005

Product Name	Anti-Human CD13 Monoclonal Antibody, Azide Free Clone B-F10
Catalog Number	854.100.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	KG-1
Specificity	Recognises aminopeptidase N, a 150-170 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells
Immunisation	KG-1