

Product Datasheet

Anti-Clostridium difficile (toxin A) V5148

Product Name	Anti-Clostridium difficile (toxin A)
Catalog Number	V5148
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with toxin B
Specificity	Toxin A of Clostridium difficile
Matched Pairs	ELISA pairs (capture-conjugate) : V5154-V5148, V5157-V5148
Shelf Life	4 years