

Product Datasheet

Anti-Dengue virus (NS1) V5162

Product Name	Anti-Dengue virus (NS1)
Catalog Number	V5162
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, LF, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	NS1 protein of all 4 serotypes
Matched Pairs	ELISA pairs (capture-conjugate) : V5162-V5163. V5162-V5164, V5162-V5165. LF pairs (capture-conjugate) : V5166-V5162
Shelf Life	4 years