

Product Datasheet

Anti-M. tuberculosis (CFP 10) V5172

Product Name	Anti-M. tuberculosis (CFP 10)
Catalog Number	V5172
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Recognizes the CFP 10 protein of M. tuberculosis
Shelf Life	Four years from date of manufacture