

# Product Datasheet

## Anti-Parvovirus B19 V5191

|                     |   |
|---------------------|---|
| Product Name        | Anti-Parvovirus B19   |
| Catalog Number      | V5191   |
| Tested Applications | ELISA, IFA  |
| Brand               | ViroStat  |
| Form/Appearance     | Purified from mouse ascites by protein A chromatography                     |
| Concentration       | ≥ 1 mg/ml   |
| Storage             | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype             | IgG1  |
| Clonality           | Monoclonal  |
| Purity              | ≥ 90%   |
| Buffer              | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide       |
| Cross Reactivity    | Non-reactive with uninfected human cells                                    |
| Specificity         | VP2 of human Parvovirus (B19)   |
| Shelf Life          | Four years from date of manufacture   |