

Product Datasheet

Anti-Clostridium difficile (GDH) V5196

Product Name	Anti-Clostridium difficile (GDH)
Catalog Number	V5196
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with mammalian cells
Specificity	Glutamate dehydrogenase of C. difficile
Matched Pairs	ELISA pairs (capture-conjugate) : V5192-V5196, V5193-V5196, V5194-V5196, V5195-V5196
Shelf Life	4 years