

Product Datasheet

Anti-h Papilloma virus L1 (HPV 16) V5301

Product Name	Anti-h Papilloma virus L1 (HPV 16)
Catalog Number	V5301
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	IHC, ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with the L1 of type 16
Specificity	L1 of type 16
Shelf Life	4 years