

# Product Datasheet

## Anti-h Papilloma virus L1 (HPV 16 & 18) V5305

Product Name	Anti-h Papilloma virus L1 (HPV 16 & 18)
Catalog Number	V5305
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with the L1 of types 16 & 18
Specificity	L1 of types 16 & 18
Shelf Life	Four years from date of manufacture