

Product Datasheet

Anti-h Papilloma virus L1 (HPV 6, 11, 16)
V5306

Product Name	Anti-h Papilloma virus L1 (HPV 6, 11, 16)
Catalog Number	V5306
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with the L1 of types 6, 11, & 16
Specificity	L1 of types 6, 11, & 16
Shelf Life	Four years from date of manufacture