

Product Datasheet

Anti-h Papilloma virus L1 V5307

Product Name	Anti-h Papilloma virus L1
Catalog Number	V5307
Tested Applications	ELISA, WB, IHC
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with the L1 of types 6, 11, 16, 18, 31, 33, 45, 52, & 58
Specificity	L1 of types 6, 11, 16, 18, 31, 33, 45, 52, & 58
Shelf Life	Four years from date of manufacture