

## Product Datasheet

### Anti-h Papilloma virus L1 V5308

<b>Product Name</b>	Anti-h Papilloma virus L1
<b>Catalog Number</b>	V5308
<b>Description</b>	Purified from mouse ascites by protein A chromatography
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	WB, IHC, ELISA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from mouse ascites by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90 %
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Reactive with the L1 of types 6, 11, 16, 18, 31, 33, 45, 52, & 58
<b>Specificity</b>	L1 of types 6, 11, 16, 18, 31, 33, 45, 52, & 58
<b>Shelf Life</b>	4 years