

Product Datasheet

Anti-Mycoplasma genitalium V5602

Product Name	Anti-Mycoplasma genitalium
Catalog Number	V5602
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Mycoplasma hominis, Mycoplasma pneumoniae, Ureaplasma urealyticum, Herpes virus type 2, Chlamydia trachomatis, Nesseria gonorrhoeae, and human cells
Specificity	M. genitalium
Matched Pairs	ELISA pairs (capture-conjugate) : V5601-V5602, V5602-V5601, V5603-V5602, V5602-V5604
Shelf Life	Four years from date of manufacture