

# Product Datasheet

## Anti-Mumps virus V6007

Product Name	Anti-Mumps virus
Catalog Number	V6007
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Measles virus, Respiratory Syncytial virus, & Parainfluenza 1, 2, or 3 viruses
Specificity	Nucleoprotein of Mumps virus
Shelf Life	Four years from date of manufacture