

# Product Datasheet

## Anti-Measles virus V6017

|                     |   |
|---------------------|---|
| Product Name        | Anti-Measles virus  |
| Catalog Number      | V6017   |
| Tested Applications | ELISA, immunofluorescence   |
| Brand               | ViroStat  |
| Form/Appearance     | Purified from mouse ascites by protein A chromatography                     |
| Concentration       | ≥ 1 mg/ml   |
| Storage             | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype             | IgG2b   |
| Clonality           | Monoclonal  |
| Purity              | ≥ 90%   |
| Buffer              | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide       |
| Cross Reactivity    | Non-reactive with Mumps virus, RSV, Para 1, Para 2, & Para 3                |
| Specificity         | Measles virus   |
| Shelf Life          | Four years from date of manufacture   |