

Product Datasheet

Anti-Legionella pneumophila V6029

Product Name	Anti-Legionella pneumophila
Catalog Number	V6029
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG3
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with LPS of all pneumophila strains tested (11)
Specificity	LPS of all pneumophila strains tested (11)
Matched Pairs	ELISA pairs (capture-conjugate) : V6029-V6029, V6027-V6029
Shelf Life	Four years from date of manufacture