

Product Datasheet

Anti-Human CD16 Monoclonal Antibody, Azide Free Clone B-E16 854.130.005

Product Name	Anti-Human CD16 Monoclonal Antibody, Azide Free Clone B-E16
Catalog Number	854.130.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Immunogen	Natural killer cells
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the low affinity Fc receptor for IgG (Fcγ RIII), a 50-65 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells