

Product Datasheet

Anti-Giardia lamblia V6291

Product Name	Anti-Giardia lamblia
Catalog Number	V6291
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	65 kDa Giardia antigen
Matched Pairs	ELISA pairs (capture-conjugate) : V6291-V6292, V6291-V6241, V6291-V6231, V6291-V6261
Shelf Life	4 years