

Product Datasheet

Anti-Salmonella V6347

Product Name	Anti-Salmonella
Catalog Number	V6347
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, LF, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Salmonella paratyphi A LPS
Matched Pairs	Self-pairs in LF immunoassay
Shelf Life	4 years