

## Product Datasheet

### Anti-Salmonella LPS core V6381

<b>Product Name</b>	Anti-Salmonella LPS core
<b>Catalog Number</b>	V6381
<b>Description</b>	Purified from mouse ascites by protein A chromatography
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, IFA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from mouse ascites by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90 %
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Reactive with all salmonella O-serogroups tested: A, B, C1, C2, D, E1, E3, E4, F, G1, G2. Non-reactive with E. coli (mix), Klebsiella, Citrobacter, Pseudomonas, Yersinia, Shigella, Proteus & Legionella
<b>Specificity</b>	Common LPS core of all Salmonellae tested
<b>Shelf Life</b>	4 years