

Product Datasheet

Anti-Candida albicans V6402

Product Name	Anti-Candida albicans
Catalog Number	V6402
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgM
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with all C. albicans strains tested (9) + C. glabrata & C. tropicales. Non-reactive with C. parapsilosis and C. krusei
Specificity	All C. albicans strains tested (9) + C. glabrata & C. tropicales
Shelf Life	4 years