

Product Datasheet

Anti-Enolase (fungal) V6492

Product Name	Anti-Enolase (fungal)
Catalog Number	V6492
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with enolase from <i>S. cerevisiae</i> and candida species: <i>C. albicans</i> , <i>C. parapsilosis</i> , <i>C. tropicales</i> , & <i>C. glabrata</i> . Other fungi not tested
Specificity	Enolase from several candida species & <i>S. cerevisiae</i>
Matched Pairs	ELISA pairs (capture-conjugate): V6492-V6491
Shelf Life	4 years