

## Product Datasheet

### Anti-Ureaplasma V6522

**Product Name** Anti-Ureaplasma

**Catalog Number** V6522

**Tested Applications** ELISA, WB, IFA

**Brand** ViroStat

**Form/Appearance** Purified from mouse ascites by protein A chromatography

**Concentration**  $\geq 1$  mg/ml

**Storage** 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.

**Isotype** IgG1

**Clonality** Monoclonal

**Purity**  $\geq 90\%$

**Buffer** 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide

**Cross Reactivity** Reactive with serovars 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12 & 14.  
Non-reactive with Mycoplasma hominis and Mycoplasma genitalium

**Specificity** Multi-banded antigen of U. parvum and U. urealyticum serovars

**Shelf Life** Four years from date of manufacture