

# Product Datasheet

## Anti-Ureaplasma V6522

Product Name	Anti-Ureaplasma
Catalog Number	V6522
Tested Applications	ELISA, WB, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with serovars 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12 & 14. Non-reactive with Mycoplasma hominis and Mycoplasma genitalium
Specificity	Multi-banded antigen of U. parvum and U. urealyticum serovars
Shelf Life	Four years from date of manufacture