

Product Datasheet

Anti-Ureaplasma V6523

Product Name	Anti-Ureaplasma
Catalog Number	V6523
Tested Applications	ELISA, WB, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with serovars 1, 3, 6, & 14 (+ 7 & 11 weakly). Non-reactive with Mycoplasma hominis and Mycoplasma genitalium
Specificity	Multi-banded antigen of several Ureaplasma species
Shelf Life	Four years from date of manufacture