

Product Datasheet

Anti-Ureaplasma V6525

Product Name	Anti-Ureaplasma
Catalog Number	V6525
Tested Applications	ELISA, WB, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with serovars 3 & 14. Non-reactive with Mycoplasma hominis and Mycoplasma genitalium
Specificity	Multi-banded antigen of several Ureaplasma serovars
Shelf Life	Four years from date of manufacture