

Product Datasheet

Anti-Mycoplasma pneumoniae V6561

Product Name	Anti-Mycoplasma pneumoniae
Catalog Number	V6561
Tested Applications	ELISA, WB, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with multiple WB bands of Mycoplasma pneumoniae. Non-reactive with uninfected cells
Specificity	Multiple WB bands of Mycoplasma pneumoniae
Shelf Life	Four years from date of manufacture