

Product Datasheet

Anti-Helicobacter pylori (flagellin) V6641

Product Name	Anti-Helicobacter pylori (flagellin)
Catalog Number	V6641
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	WB, ELISA, LF, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive in ELISA with C. jejuni, E. coli (mix), Salmonella (mix), Shigella (mix), P. aeruginosa, Yersinia, and Citrobacter
Specificity	54 KDa (flagellin) of H. pylori
Matched Pairs	ELISA pairs (capture-conjugate): V6641-V6642, V6641-V6643, V6641-V6644
Shelf Life	4 years