

Product Datasheet

Anti-Helicobacter pylori (HSP) V6661

Product Name	Anti-Helicobacter pylori (HSP)
Catalog Number	V6661
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	WB, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with the 58 kDa (HSP) of <i>H. pylori</i> . Non-reactive in ELISA with <i>C. jejuni</i> , <i>E. coli</i> (mix), <i>Salmonella</i> (mix), <i>Shigella</i> (mix), <i>P. aeruginosa</i> , <i>Yersinia</i> , and <i>Citrobacter</i>
Specificity	58 kDa (HSP) of <i>H. pylori</i>
Shelf Life	4 years