

Product Datasheet

Anti-Campylobacter jejuni V7721

Product Name	Anti-Campylobacter jejuni
Catalog Number	V7721
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	WB, ELISA, LF
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with a ~ 100kDa antigen in WB. Non-reactive with H. pylori, E. coli (several), Salmonellae (several), Shigellae (several), P. aeruginosa, Yersinia & Citrobacter
Specificity	C. jejuni
Shelf Life	4 years