

Product Datasheet

Anti-Staphylococcus Enterotoxin (toxin B) V7931

| | |
|----------------------------|---|
| Product Name | Anti-Staphylococcus Enterotoxin (toxin B) |
| Catalog Number | V7931 |
| Tested Applications | ELISA |
| Brand | ViroStat |
| Form/Appearance | Purified from mouse ascites by protein A chromatography |
| Concentration | ≥ 1 mg/ml |
| Storage | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥ 90% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Not determined |
| Specificity | Enterotoxin B of <i>S. aureus</i> |
| Matched Pairs | ELISA pairs (capture-reporter) V7941-V7931 |
| Shelf Life | Four years from date of manufacture |