

# Product Datasheet

## Anti-Staphylococcus Enterotoxin (toxin B) V7941

Product Name	Anti-Staphylococcus Enterotoxin (toxin B)
Catalog Number	V7941
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Enterotoxin B of S. aureus
Matched Pairs	ELISA pairs (capture-reporter) V7941-V7931
Shelf Life	Four years from date of manufacture