

# Product Datasheet

## Anti-Vaccinia virus V8116

Product Name	Anti-Vaccinia virus
Catalog Number	V8116
Tested Applications	ELISA, IHC, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with free virus, uninfected cells, Lister and NYCBH strains of vaccinia
Specificity	A33R protein
Shelf Life	Four years from date of manufacture