

Product Datasheet

Anti-Bovine Papillomavirus V5091

Product Name	Anti-Bovine Papillomavirus
Catalog Number	V5091
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with HPV16-L1, HPV11-L1 and uninfected bovine cells
Specificity	Bovine papillomavirus
Shelf Life	Four years from date of manufacture