

Product Datasheet

Anti-Feline Leukemia virus (p27) V7218

Product Name	Anti-Feline Leukemia virus (p27)
Catalog Number	V7218
Description	Purified from mouse ascites by protein A chromatography
Species/Host	Mouse
Tested Applications	WB, LF, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	FelV p27
Matched Pairs	Pairs for p27 detection (capture-conjugate) : V7226-V7218, V7228-V7218
Shelf Life	4 years