

Product Datasheet

Anti-Canine Distemper virus V7516

Product Name	Anti-Canine Distemper virus
Catalog Number	V7516
Description	Purified from mouse ascites by protein A chromatography
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	+2-8 °C or -20 °C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Canine Parvovirus and Canine Coronavirus
Specificity	Nucleocapsid of CDV
Matched Pairs	ELISA pairs (capture-conjugate) : V7516-V7515
Shelf Life	4 years