

Product Datasheet

Dengue virus (Env, Serotypes 1-4) Antigen, Recombinant V8824

Product Name	Dengue virus (Env, Serotypes 1-4) Antigen, Recombinant
Catalog Number	V8824
Species/Host	E. coli
Tested Applications	ELISA, LF
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	6X-His
Purity	~95% by SDS-PAGE
Buffer	PBS, pH 7.4
Specificity	Envelope antigen, contains serotype common regions
Molecular Weight	19 kDa