
Product Datasheet

Zika virus (Env) Antigen, Recombinant V8826

Product Name	Zika virus (Env) Antigen, Recombinant
Catalog Number	V8826
Species/Host	E. coli
Tested Applications	ELISA, LF
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	6X-His - C terminus
Purity	~95% by SDS-PAGE
Buffer	PBS + arginine + 0.01% azide
Specificity	Envelope antigen, truncated