
Product Datasheet

Zika virus (Env & NS1) Antigen, Recombinant V8827

Product Name Zika virus (Env & NS1) Antigen, Recombinant

Catalog Number V8827

Species/Host E. coli

Tested Applications ELISA, LF

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note 6X-His tag

Purity ~90% by SDS-PAGE

Buffer 20 mM Phosphate, 8.0 + 10 mM NaCl

Specificity Envelope-NS1 chimera, proprietary

Molecular Weight 15 kDa