

Product Datasheet

Hepatitis C virus (core) Antigen, Recombinant V8903

Product Name	Hepatitis C virus (core) Antigen, Recombinant
Catalog Number	V8903
Species/Host	E. coli
Tested Applications	WB, ELISA, LF
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	β-galactosidase (114 kDa)
Purity	~95% by SDS-PAGE
Buffer	0.05M Tris-saline, pH 8.0 + 6M urea
Specificity	Core antigen of HCV, truncated from C-terminus
Molecular Weight	133 kDa