

## Product Datasheet

### Human Immunodeficiency Virus 1 (p24) Antigen, Recombinant V8909

<b>Product Name</b>	Human Immunodeficiency Virus 1 (p24) Antigen, Recombinant
<b>Catalog Number</b>	V8909
<b>Species/Host</b>	Pichia pastoris
<b>Tested Applications</b>	WB, ELISA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified recombinant antigen
<b>Storage</b>	≤ -20°C
<b>Note</b>	β-galactosidase (114 kDa)
<b>Purity</b>	~95%
<b>Buffer</b>	0.02 M Sodium phosphate, pH 7.2
<b>Specificity</b>	Gag, p24
<b>Molecular Weight</b>	148 kDa