

---

## Product Datasheet

### Hepatitis B virus Core (HBcAg) Antigen, Recombinant V8914

<b>Product Name</b>	Hepatitis B virus Core (HBcAg) Antigen, Recombinant
<b>Catalog Number</b>	V8914
<b>Species/Host</b>	E. coli
<b>Tested Applications</b>	WB, ELISA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified recombinant antigen
<b>Storage</b>	≤ -20°C
<b>Note</b>	No fusion partner present
<b>Purity</b>	~95% by SDS-PAGE
<b>Buffer</b>	50 mM Tris, 1mM EDTA, pH 7.5 + 150 mM NaCl
<b>Specificity</b>	Core antigen, 180 a.a.
<b>Molecular Weight</b>	17 kDa