

Product Datasheet

Human Immunodeficiency Virus 1 (nef) Antigen, Recombinant V8941

Product Name	Human Immunodeficiency Virus 1 (nef) Antigen, Recombinant
Catalog Number	V8941
Species/Host	E. coli
Tested Applications	WB, ELISA
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	β-galactosidase
Purity	~95% by SDS-PAGE
Buffer	10 mM Carbonate + 10 mM BME + 10 mM EDTA
Specificity	Nef, a.a.3 to a.a. 190