

Product Datasheet

Treponema pallidum (TmpA) Antigen, Recombinant V8948

Product Name	Treponema pallidum (TmpA) Antigen, Recombinant
Catalog Number	V8948
Species/Host	Pichia pastoris
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	No fusion partner present
Purity	~95% by SDS-PAGE
Buffer	20 mM Tris, pH 8.0 + 8M urea + 10 mM BME
Specificity	TmpA, full length
Molecular Weight	14 kDa