

Product Datasheet

Anti-Human CD21 Monoclonal Antibody, FITC Conjugated Clone B-E5 854.161.010

Product Name	Anti-Human CD21 Monoclonal Antibody, FITC Conjugated Clone B-E5
Catalog Number	854.161.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Daudi cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 5% BSA and 0.1 % sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Epstein-Barr virus receptor (EBV R), CR2, C3d R, a 110 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells