

Product Datasheet

Feline Immunodeficiency Virus (core) Antigen, Recombinant V8963

Product Name	Feline Immunodeficiency Virus (core) Antigen, Recombinant
Catalog Number	V8963
Species/Host	E. coli
Tested Applications	ELISA, LF
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	2-8°C
Note	Maltose binding protein (truncated)
Purity	~90% by SDS-PAGE
Buffer	Bicarbonate buffer, pH 8.5 + 0.02% azide
Specificity	Core
Molecular Weight	37 kDa