

Product Datasheet

Plasmodium vivax (MSP) Antigen, Recombinant V8969

Product Name	Plasmodium vivax (MSP) Antigen, Recombinant
Catalog Number	V8969
Species/Host	E. coli
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	Maltose binding protein (truncated)
Purity	~95% by SDS-PAGE
Buffer	PBS, pH 7.2 + 0.1% sodium azide
Specificity	MSP protein, 108 a.a.
Molecular Weight	~ 66 kDa