

Product Datasheet

Hepatitis C virus (NS3) Antigen, Recombinant V8995

Product Name	Hepatitis C virus (NS3) Antigen, Recombinant
Catalog Number	V8995
Species/Host	E. coli
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	2-8°C
Note	His tag
Purity	~95% by SDS-PAGE
Buffer	50 mM Tris saline, pH 8.5
Specificity	NS3 antigen of HCV, C-terminal 380 aa
Molecular Weight	16 kDa