

## Product Datasheet

### Human Immunodeficiency Virus 1 (Subtype B gp 41) Antigen, Recombinant V9000

<b>Product Name</b>	Human Immunodeficiency Virus 1 (Subtype B gp 41) Antigen, Recombinant
<b>Catalog Number</b>	V9000
<b>Species/Host</b>	E. coli
<b>Tested Applications</b>	LF, ELISA, WB
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified recombinant antigen
<b>Storage</b>	≤ -20°C
<b>Note</b>	His tag
<b>Purity</b>	~95% by SDS-PAGE
<b>Buffer</b>	20 mM Na-carbonate, pH 9.6
<b>Specificity</b>	Ectodomain Subtype B, gp41
<b>Molecular Weight</b>	13 kDa